



TechLITE Pen Needle Compatibility Chart

TechLITE Pen Needles fit most standard insulin pens.

Standard Needle-based Injection System EN ISO 11608-2:2012 clause 4.9 and 11	
Manufacturer / Distributor:	Brand Name:
B. Braun	Omnican®Pen 31
	Omnican®Pen 32
BD Medical	PuregonPen®
Copernicus Sp. z o.o. / Bioton	GensuPen®
Diamesco Co., Ltd./pendiq GmbH	pendiq 2.0
Gani Technology	Xiu Lin Pen
Eli Lilly	Humalog®KwikPen 100 U/ml, 200 U/ml
	HumaPen® Ergo
	HumaPen® Ergo II
	HumaPen® L U X U R A
	HumaPen® L U X U R A HD
	HumaPen® MEMOIR
	HumaPen® savvio
	Humulin® N Pen 100U/ml
	βyetta® Pen
	ABASAGLAR® KwikPen® 100 Units/ml
	Saline KwikPen™
	Forsteo® 20 µg/80µl
	Haselmeier
i-penⓂ	
Haselmeier / Berlin-Chemie	BerliPen® JUNIOR
	BerliPen® areo2
	BerliPen® PRECISION
	BerliPen® areo
	JuniorSTAR®
	TactiPen®

Standard Needle-based Injection System EN ISO 11608-2:2012 clause 4.9 and 11	
Manufacturer / Distributor:	Brand Name:
Novo Nordisk	FlexPen®
	FlexTouch®
	NovoLet®
	NovoPen® 3
	NovoPen® 4
	NovoPen® 5
	NovoPen Echo®
	NovoPen® Junior
	Xultophy®
	Victoza® 6mg/ml
Owen Mumford	Tresiba® FlexTouch® 200 Units/mL
	Levemir®InnoLet® 100 Units/ml
	Autopen® 24 3.0ml x 1 unit 1-21units
	Autopen® 24 3.0ml x2 units 2-42units
	Autopen® Classic 3.0ml x 1 unit 1-21units
Owen Mumford / Berlin-Chemie	Autopen® Classic 3.0ml x2 units 2-42units
	BerliPen® 301
PENDIQ	BerliPen® 302
	pendiq
Sanofi-Aventis	ClikSTAR®
	OptiPen® Pro1
	OptiPen® Pro2
	SoloSTAR®
	Lantus® SoloStar® 100 Units/ml
	Lyxumia® 20 µg
Ypsomed	ypsopen®
Ypsomed / Фармстандарт	БиоматикПен®
Ypsomed / Tonghua Dongbao	Gansulin Pen

NOTE! The list above is not exhaustive - the brand names may vary on different markets, change due to mergers of acquisitions of their owners and private label arrangements. All specifications are subject to change. Consult appropriate company for additional information. Trademarks are the property of their respective owners.